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- ¹ Qualification of High-Recovery, Flocked Swabs as Compared to Traditional Rayon Swabs for Microbiological Environmental Monitoring of Surfaces – Giblerto Dalmaso, Manuela Bini, Roberto Paroni and Michela Ferrari
- ² New Sampling Approach to Detect Environmental Microbial Contamination in Italian Hospital Wards - Marco Ferrari, Adriano Anesi Unit of Hospital Hygiene, Laboratory of Microbiology, ASST of Lodi, Italy
- ³ Surveillance Program of Environmental Microbial Contamination in Italian Health Care Facilities – Marco Ferrari, Matteo Beffa, Adriano Anesi Unit of Hospital Hygiene, University of Brescia, Italy, Laboratory of Microbiology, ASST of Lodi Italy
- ⁴ Strengthening surveillance of and response to foodborne diseases: a practical manual. Investigating foodborne disease outbreaks
- ⁵ The Review on Antimicrobial Resistance – The global challenge of drug-resistant infections – Hala Audi
- ⁶ Changing epidemiology of *Listeria monocytogenes* outbreaks, sporadic cases, and recalls globally: A review of ProMED reports from 1996 to 2018
- ⁷ Only for section regarding swab stick treatment/sampling/transport

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> European Patent #1608268
 > US Patent #8,114,027
 > US Patent #8,317,728
 > US Patent #8,979,784

> US Patent #9,011,358
 > US Patent #9,173,779
 > Canadian Patent #2,515,205
 > Japanese Patent # 2007-523663

> Australian Patent #2004226798
 > New Zealand Patent #541560
 > Chinese Patent #ZL200610099310.9

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APPLIED AND ENVIRONMENTAL MICROBIOLOGY, Aug. 2010, p. 5148–5158.
0099-2240/10/512.00 doi:10.1128/AEM.00399-10
Copyright © 2010, American Society for Microbiology. All Rights Reserved.

Vol. 76, No. 15

Validation of a Nylon-Flocked-Swab Protocol for Efficient Recovery of Bacterial Spores from Smooth and Rough Surfaces

Alexander Probst,¹ Rainer Facius,² Reinhard Wirth,¹ and Christine Moissl-Eichinger^{1*}

¹Institute of Microbiology and Archaea Center, University of Regensburg, Universitaetsstrasse 31, 93053 Regensburg, Germany and
²German Aerospace Center (DLR), Linder Hoehe, 51147 Cologne, Germany

Received 12 February 2010/Accepted 1 June 2010

In order to meet planetary-protection requirements, culturable bacterial spore loads are measured representatively for the total microbial contamination of spacecraft. However, the National Aeronautics and Space Administration's (NASA's) cotton swab protocols for spore load determination have not changed for decades. To determine whether a more efficient alternative was available, a novel swab was evaluated for recovery of different *Bacillus atrophaeus* spore concentrations on stainless steel and other surfaces. Two protocols for the nylon-flocked swab (NFS) were validated and compared to the present NASA standard protocol. The results indicate that the novel swab protocols recover 3- to 4-fold more (45.4% and 49.0% recovery efficiency) *B. atrophaeus* spores than the NASA standard method (13.2%). Moreover, the nylon-flocked swab protocols were superior in recovery efficiency for spores of seven different *Bacillus* species, including *Bacillus anthracis* Sterne (recovery efficiency, 20%). The recovery efficiencies for *B. atrophaeus* spores from different surfaces showed a variation from 5.9 to 62.0%, depending on the roughness of the surface analyzed. Direct inoculation of the swab resulted in a recovery rate of about 80%, consistent with the results of scanning electron micrographs that allowed detailed comparisons of the two swab types. The results of this investigation will significantly contribute to the cleanliness control of future life detection missions and will provide significant improvement in detection of *B. anthracis* contamination for law enforcement and security efforts.

The recent discovery of liquid water on Mars has sparked debate about the possibility of extraterrestrial life (37). Consequently, highly sensitive biosensors will be deployed onboard spacecraft like the Mars Science Laboratory (MSL), using technologies such as gas chromatographical analysis to search for the smallest traces of life (<http://mars.jpl.nasa.gov/msl/mission/>). Contamination of equipment by terrestrial microorganisms resulting from a lack of spacecraft cleanliness could significantly compromise the integrity of life detection missions and result in falsely positive extraterrestrial life signals. The prevention of this so-called "forward contamination" is one major goal of American and European space agencies' planetary-protection efforts. Regular determination of a spacecraft's bioload and the mission components throughout assembly are mandatory for detecting unacceptably high contamination that exceeds levels set by the United Nations treaty (Outer Space Treaty [11]).

Modern spacecraft hardware is very susceptible to standard heat sterilization protocols, so baking the entire spacecraft, such as the Viking Lander Capsule at 111.7°C ± 1.7°C for 23 to 30 h is no longer feasible (30). Alternative cleaning and sterilization methodologies for spacecraft components prior to assembly (i.e., nonthermal plasma technologies) have been discussed (36). However, after integration, sterile hardware is

limited access to microbial cleanliness of

Institute of Microbiology
Universitaetsstrasse 31,
93053 Regensburg, Germany
Tel: 49 941 943-3100
Fax: 49 941 943-3101
Email: am-i-regensburg.de

5148

a spacecraft and its surroundings is meticulously maintained through frequent cleaning and sterilization routines. The regular and frequent detection of possible contamination in the assembly environment is more important than ever.

To estimate the severity of microbial contamination, the National Aeronautics and Space Administration's (NASA) standard procedure focuses on aerobic spore-forming bacteria. Briefly, surface samples are taken from spacecraft components with cotton swabs or wipes. After an extraction procedure, the samples are subjected to a short heat shock (15 min at 60°C) to kill vegetative cells and then pour plated in *Trypsin Soy* (TSA) for the enumeration of CFU. This protocol was originally developed for the Viking mission more than 30 years ago (30) and has remained, for the most part, unchanged.

Recent studies have shown that cotton swabs have acceptable recovery efficiencies for *Bacillus* spores (45.4% and 49.0%) due to their organic nature, may leave residues on smooth surfaces. Furthermore, their comparatively high moisture content could lead to false positives or inhibition shown in NASA one day incorporate molecular technologies into their microbial-detection protocols (7).

Based on these observations, researchers are beginning to move away from cotton in favor of alternative swabs made from rayon or macrofoam (6, 18). A recent study reported high recovery efficiencies for various vegetative cell types on stainless steel surfaces by applying a novel swab with a bulb-shaped head flocked with nylon fibers (12). Patented in 2004, this design facilitates the release of particles and microbes, resulting in a significantly higher detection rate. The broad applicability of these nylon-flocked swabs (NFS) has been demonstrated by their use in various situations, including isolating pathogens from medical environments (10, 20).

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